

Doffing Operating Room PPE

STEP 1

Grab yellow isolation gown at waist and pull forward to break ties at back. Continue pulling forward to cause back of gown to tear away.

STEP 2

Roll gown into itself and peel off gloves at same time. Only touch inside of gown and gloves. Discard and perform hand hygiene.

STEP 3

Remove face shield or eye goggles from back without touching face. Consider saving for cleaning or re-use. Perform hand hygiene.

STEP 4

Remove N-95 by lifting bottom strap and then top strap over head. If not visibly soiled, consider continuous use or re-use.

USE "BUDDY SYSTEM" TO MONITOR FOR PROPER TECHNIQUE

CHANGE YOUR SCRUBS PRIOR TO PROVIDING FURTHER CLINICAL CARE OR LEAVING HOSPITAL