

Cardiac Room Setup

On Mayo Stand

- PSI Kit
- Mayo stand sterile cover



In Omnicell Drug Drawer

- Phenylephrine infusion bag – 80mcg/mL (20mg in 250mL of NS)
- Nitroglycerin push - 40mcg/mL in 10mL syringe (1mL of 400mcg/mL solution diluted in 9mL of NS)
- Esmolol push - 10mg/mL in 10mL syringe
- Norepinephrine push – 1.6mcg/mL (take 4mg of norepinephrine and place in 250mL NS bag; then dilute 1mL of that in 9mL of NS)
- Premade phenylephrine and ephedrine bolus syringes
- Premade norepinephrine syringe - 80mcg/mL



On Anesthesia Tray

- Mac 3/Miller 2 blades with handle
- 7.5 ETT with 10mL syringe (unopened)
- #3 and #5 oral airways
- Temp probe
- Nasogastric tube
- Benzoin bullet x 2
- Lubricating jelly x 2
- Alcohol swabs
- SedLine EEG monitor
- EKG red dots and pulse oximeter sticker
- Anesthesia Machine check
- Working suction



On Omnicell

- Arterial line kit (arm board with taped OB sponges, sterile 4x4, cloth and plastic tape, prepared)



- flush, lidocaine syringe, benzoin bullet, 20G angiocath x 2
- IV kit with 14 & 16G angiocath and lidocaine syringe
- Nitroglycerine bottle (250mL)

On OR Table

- Standard monitor cables and arterial, CVL, PAC transducer setup (techs set this up, but ensure it's set up)

On IV Pole

- 1L LR bag on pump set
- 500mL NS bag on buretrol

Next to Omnicell (resting on sharps container)

- Transesophageal echo probe

